

Remember, With Home Buddy You're Never Home Alone!

Medication Management

How it works...

Alarm sounds when dose is ready, if not taken, unit will contact Care Center.

1



Simple & Reliable

1. Load
2. Set clock
3. Program alarm
4. Close

Here's how it works...

When the time comes to take your next dose, the alarm will sound and the dose will rotate into position. The alarm stops when the dose is taken out.

If the dose is not removed within 30 minutes, the unit can contact the 24-hour Care Center which then attempts to alert someone of a non-compliance event.

It's really that easy!

ALWAYS
24/7
ACTIVE



No hard to open pill bottles or boxes.

Home Buddy, LLC
1965 Pleasant View Circle
Wichita, KS 67203

every dose...on time

2 “Ed this is the Care Center it's time to take your medication.”



3 “Oh, OK thanks for the reminder, I guess I didn't hear it and was so busy I didn't realize it was time to take my pills out of the machine. Thanks so much, I'll go in and take them right now.”

Lots of Prescriptions?

Keeping up with a lot of prescriptions can be a real challenge. There may be one set of pills and vitamins to be taken in the morning, another set to be taken throughout the afternoon, and still other pills before bed.

Now there's *an easier way* to manage all those pills. The Medready simplifies everyday medication management, giving you more time and freedom. It is simple and reliable, and the convenient medication dispenser is much easier to open than a pill box or prescription bottle.

Our goal has always been to keep people safe and independent as long as possible. Helping people manage their prescription therapies more effectively is key to accomplishing that goal.

**Right Med
Right Dose
Right Time**

You'll have fewer worries about...

- Forgetting to take a dose
- Taking a double dose
- Taking at the wrong time
- Identifying pills or pill bottles for people with decreased vision
- Opening and closing pill bottles
- Privacy regarding the type and quantity of medication